



Shuster Laboratories, Inc.  
85 John Road  
Canton, MA 02021  
Phone: 781-821-2200  
Fax: 781-821-9266  
www.shusterlabs.com

**TAQA®**  
**(Technically Advanced Quality Assurance)**  
**Evaluation Report**

Mr. Mahmoud M. Sales  
**SCI International, Inc.**  
5902 Enterprise Court  
Fredrick, MD 21703

Folder #: 0807895  
Date: 4/30/2008  
Received Date: 4/9/08  
Page: 1 of 3

Product:	New Choice Ovulation Cassette Tests	UPC Code:	
Size:	Sample	Reference #:	Partial Prepaid Chk #2637
Supplier:	SCI International (Made in Canada)	Inventory ID #:	143376
Code:	Lot L09256M5 EXP 01/2009	Lab Sample #:	0807895001
Purpose:	TAQA EVALUATION vs. CLEARBLUE	--	

**National Reference Brand:** Clearblue Easy Blue Ovulation Test; Lot Code: LOT NO: 7310937291 EXP DATE: 2009-07 - provided by SCI International

**Overall Rating:** 10.0 Sample provided meets TAQA quality criteria for the parameters examined. See note.

**Comments:** Note: SCI requested that this product be tested for overall performance against the ClearBlue "stick" product as there is no longer a National Brand cassette type product available at the current time. SCI is anticipating that this type of product will be in demand in the future.

Overall Quality Rating \*

- 9.0 - 10.0 Sample tested is Comparable to the National Brand; meets TAQA quality criteria
- 7.5 - < 9.0 Sample tested is Marginally Comparable to the National Brand; does not meet all TAQA quality criteria
- < 7.5 Sample tested is Not Comparable to the National Brand; does not meet TAQA quality criteria

Maryann W. Sullivan  
Client Services Manager

Determination		Clearblue Easy Blue Ovulation Test	Sample of SCI New Choice Ovulation Cassette Tests (SCI)	Comments
<b>Performance – Bench Top</b>	Positive Result	Obvious positive result 38 seconds after 20 second dip; result remained after 10 minutes; meets claim of within 3 minutes	Obvious positive result 53 seconds after 20 second dip; result remained after 10 minutes; meets claim of within 10 minutes	
	Negative Result	Obvious negative result 60 seconds after 20 second dip; result remained after 10 minutes; meets claim of within 3 minutes	Obvious negative result 50 seconds after 20 second dip; result remained after 10 minutes; meets claim of within 10 minutes	
<b>Physical</b>	Cassette Dimensions, inches	--	--	
	Width	--	11/16	
	Thickness	--	0.22	
	Overall Length	--	3-3/16	
	Product Description	Sticks are individually wrapped in a heat-sealed foil package; instruction pamphlet included; boxed	Rectangular cassette containing 2 recessed windows: a rectangular one to display results, a circular one for urine drop from plastic dropper; both are packaged in a heat-sealed foil package	Opened carton included for use instructions

		Clearblue Easy Blue Ovulation Test	Sample of SCI New Choice Ovulation Cassette Tests	Comments
<b>Physical</b>	Sample Well	--	Circular recessed window	
	Windows	--	One recessed rectangular window; one recessed circular window	
	Ease of Opening/ Closing Container	--	Acceptable	
	Tamper Evident Feature	--	Sample packaging	
	Product Workmanship	--	Acceptable	
	Package Workmanship	--	Sample packaging	
	Package Description	--	Sample packaging	