



Specialized Technology Resources
 85 John Road
 Canton, MA 02021
 Phone: 781-821-2200
 Fax: 781-821-9266
 www.STRQuality.com

TAQA®
Technically Advanced Quality Assurance TAQA®
Shuster Report

Mr. Mahmoud M. Sales
SCI International, Inc.
 5902 Enterprise Court
 Fredrick, MD 21703

Folder #: 0912720
 Date: 5/29/2009
 Received Date: 4/7/09
 Page: 1 of 3

cc:

Product:	SCI Brand Fertility Test - Midstream*	UPC Code:	
Size:	Sample	Reference #:	
Supplier:	SCI International, Inc.	Inventory ID #:	166202
Code:	LOT L04098M11 EXP 04/2010	Lab Sample #:	0824663001
Purpose:	TAQA EVALUATION vs. FIRST RESPONSE FERTILITY TEST		

National Reference Brand: First Response Fertility Test for Women; Lot Code: BU8357FA EXP 12-2010

Overall Rating: 10.0 Sample submitted meets TAQA quality criteria for parameters examined

Note: * SCI Brands includes: New Choice, Perfect Choice, Baby Test, Baby Check, UCheck

Overall Quality Rating *

- 9.0 - 10.0 Sample tested is Comparable to the National Brand; meets TAQA quality criteria
- 7.5 - < 9.0 Sample tested is Marginally Comparable to the National Brand; does not meet all TAQA quality criteria
- < 7.5 Sample tested is Not Comparable to the National Brand; does not meet TAQA quality criteria

Maryann W. Sullivan

Maryann W. Sullivan, Client Services Manager

Determination		First Response Fertility Test for Women	SCI Fertility Test - Midstream (See note on cover)	Comments
Performance - Bench Top	Positive Results Used FSH (+) Control 50miu/mL 0912720 185488 Exp: 5/27/09	Using a 5 second dip "urine" entered results window at 14 seconds; first reaction at 70 seconds; complete results at 85 seconds; constant for 10+ min	Using a 5 second dip "urine" entered results window at 12seconds; first reaction at 3:40 minutes; complete results at 5 minutes; constant for 10+ min	
	Negative Results Used FSH (-) Control 50miu/mL 0912720 185488 Exp: 5/27/09	Using a 5 second dip "urine" entered results window at 15 seconds; first reaction at 70 seconds; complete results at 100 seconds; constant for 10+ min	Using a 5 second dip "urine" entered results window at 15 seconds; first reaction at 45 seconds; complete results at 82 seconds; constant for 10+ min	
Physical	DIMENSIONS, IN. Overall Length Width at Cap Width at Handle	5-5/8 13/16 5/8	5-5/8 5/8 5/8	
	ABSORBENT TIP Length Width Thickness	1-1/16 7/16 0.10	1-3/4 3/8 1/8	
	HANDLE Length Width Thickness	4-3/8 13/16 5/16	3-1/2 5/8 5/16	
	CAP Length Width Thickness	1-5/8 1-3/16 3/8	2-1/8 5/8 5/16	

		First Response Fertility Test for Women	SCI Fertility Test - Midstream (See note on cover)	Comments
Physical	Construction	Opaque/white plastic contoured handle with clear/colorless plastic tapered cap; white rectangular-shaped absorbent tip at one end. One oval result window covered with plastic at mid-section of handle; thumb grips at end of handle; First Response and symbol embossed on lower handle.	Opaque/white plastic tapered stick with cover; stick has an extruded (plastic) sampling collector "tray" containing an absorbent strip at one end, a rectangular result view window at mid section and a depressed finger hold at the other end. (samples submitted)	
	Product Workmanship	Acceptable	Acceptable	
	Packaging	Individually wrapped in printed heat-sealed foil pouch; boxed with clear plastic overwrap. Instruction pamphlet with illustrations inside.	Individually wrapped sticks in foil packets (used First Response instructions) (samples submitted)	