



Specialized Technology Resources  
 85 John Road  
 Canton, MA 02021  
 Phone: 781-821-2200  
 Fax: 781-821-9266  
 www.STRQuality.com

## Technically Advanced Quality Assurance TAQA® Evaluation Report

Mr. Mahmoud M. Sales  
**SCI International, Inc.**  
 5902 Enterprise Court  
 Fredrick, MD 21703

Folder #: 1013383-2 & 10152170  
 Date: 6/25/2010 & 7/23/2010  
 Received Date: 6/8/10 & 7/1/2010  
 Page: 1 of 3

cc: Mr. Stuart Neuman

Product:	SCI Brand Perfect Choice Pregnancy Midstream Test*	UPC Code:
Size:	Sample	Reference #: Prepaid American Express
Supplier:	SCI International, Inc.	Inventory ID #: 239461 & 242896*
Code:	Lot 06119M11 Exp 12/2011 and Lot 12159M11 Exp 07/2012*	Lab Sample #: 1013383001 & 1015217002*
Purpose:	TAQA EVALUATION vs. E.P.T. EARLY RESPONSE	

**National Reference Brand:** E.P.T. Early Pregnancy Test Result, (±) Results 5 days Before Missed Period;  
 Code: Exp 2012 04 Lot No. 27092 - provided by SCI International, Inc.

**Overall Rating:** 10.0 Sample submitted meets TAQA quality criteria for parameters examined (a).

**Comments;** (a) Differs in construction from E.P.T. but result/performance is comparable.

**Note:** \* SCI Brands includes: New Choice, Perfect Choice, Baby Test, Baby Check, UCheck

### Overall Quality Rating \*

- 9.0 - 10.0** Sample tested is Comparable to the National Brand; meets TAQA quality criteria
- 7.5 - < 9.0** Sample tested is Marginally Comparable to the National Brand; does not meet all TAQA quality criteria
- < 7.5** Sample tested is Not Comparable to the National Brand; does not meet TAQA quality criteria

Digitally signed by Maryann  
 Sullivan  
 Date: 2010.07.23 14:15:38 -04'00'

Maryann W. Sullivan, Client Services Manager

Determination		E.P.T. Early Pregnancy Test (±) Results 5 days Before Missed Period	SCI Brand Perfect Choice Pregnancy Midstream Test	Comments
<b>Performance - Bench Top</b>  <b>Used Sample Inv. # 242895</b>	<b>Positive Results</b>  Used HCG Hormone Positive Control 20miu/mL Lot 060710P Exp 6/20/10 (supplied by SCI)	Using a 20 sec. dip “urine” entered the round window at 13 sec. with a first reaction at 28 sec.; “urine” entered square window at 29 sec. with results at 41 seconds; complete positive at 53 sec.; constant for 10 minutes	Using a 5 second dip “urine” entered the window at 12 seconds; first sign of a reaction at 28 seconds; complete positive result at 38 seconds; constant for 10 minutes	Used SCI’s instructions
	<b>Negative Results</b>  Used Quantimetrix The Dipper Lot: 44501 Exp 2011-01	Using a 20 second dip “urine” entered round window at 20sec. with first reaction at 29 sec.; “urine” entered square window at 29sec. with a results at 44 sec.; complete negative results at 55 sec.; constant for 10 minutes	Using a 5 second dip “urine” entered the window at 12 seconds; first sign of a reaction at 30 seconds; complete positive result at 40 seconds; constant for 10 minutes	Used SCI’s instructions
<b>Physical</b>	<b>DIMENSIONS, inches</b> Overall Length Width at Cap Width at Handle	5-1/2 1 7/8	5-5/8 5.8 3/8 (at tip)	
	<b>ABSORBENT TIP</b> Length Width Thickness	1 (rounded end) 1/2 1/8	1-3/4 (at tapered end of stick) 3/8 1/8	
	<b>HANDLE</b> Length Width Thickness	4-3/8 1 (widest) - 5/16 (tapered) 3/8 (widest)	3-1/2 5/8 - 3/8 (tapered) 5/16	
	<b>CAP</b> Length Width Thickness	2-1/4 - 1-3/4 (tapered) 5/8 - 7/8 slightly tapered 5/16	2-1/8 (tapered) 5/8 - 3/8 (tapered) 5/16	

		E.P.T. Early Pregnancy Test (±) Results 5 days Before Missed Period	SCI Brand Perfect Choice Pregnancy Midstream Test	Comments
<b>Physical</b>	Construction	Opaque/white plastic contoured handle with opaque/purple plastic tapered cap; white rectangular-shaped absorbent tip with rounded end. At mid section, one oval result window covered with plastic with 2 sections: one round, one rectangular; thumb grips line near tip end, several raised lines at end of handle; "e.p.t." printed in purple on handle	Opaque/white plastic tapered stick with purple tapered plastic cover; stick has an extruded (plastic) sampling collector "tray" containing an absorbent strip at one end, a rectangular result view window at mid section and a depressed finder hold at the other end. (a) (samples submitted)	(a) Differs in construction from E.P.T. but result/performance is comparable
	Product Workmanship	Acceptable	Acceptable	
	Packaging	Individually wrapped in printed heat-sealed foil pouch; boxed with clear plastic overwrap. Instruction pamphlet with illustrations inside.	Individually wrapped sticks in foil heat-sealed packets (use instructions submitted)	