



Specialized Technology Resources
 85 John Road
 Canton, MA 02021
 Phone: 781-821-2200
 Fax: 781-821-9266
 www.STRQuality.com

Technically Advanced Quality Assurance TAQA® Evaluation Report

Mr. Mahmoud M. Sales
SCI International, Inc.
 5902 Enterprise Court
 Fredrick, MD 21703

Folder #: 1013376-1 & 1015170
 Date: 6/25/2010 & 7/23/2010
 Received Date: 6/8/10 & 7/1/10
 Page: 1 of 3

cc: Mr. Stuart Neuman

Product:	SCI Brand New Choice Pregnancy Midstream Test*	UPC Code:	
Size:	Sample	Reference #:	Prepaid American Express
Supplier:	SCI International, Inc.	Inventory ID #:	239460 and 242895*
Code:	Lot 01209C03 Exp 04/2011 and Lot 06280C03 Exp 07/2012*	Lab Sample #:	1013376001 & 1015170*
Purpose:	TAQA EVALUATION vs. FIRST RESPONSE		

National Reference Brand: First Response Pregnancy Test Early Result (6 days before a missed period);
 Lot Code: BU0068PA Exp 02-2012 - provided by SCI International

Overall Rating: 10.0 Sample submitted meets TAQA quality criteria for parameters examined.

Note: * SCI Brands includes: New Choice, Perfect Choice, Baby Test, Baby Check, UCheck

Overall Quality Rating *

- 9.0 - 10.0** Sample tested is Comparable to the National Brand; meets TAQA quality criteria
- 7.5 - < 9.0** Sample tested is Marginally Comparable to the National Brand; does not meet all TAQA quality criteria
- < 7.5** Sample tested is Not Comparable to the National Brand; does not meet TAQA quality criteria

Digitally signed by Maryann Sullivan
 Date: 2010.07.23 14:04:27 -04'00'

Maryann W. Sullivan, Client Services Manager

Determination		First Response Pregnancy Test Early Result (6 days before a missed period)	SCI Brand New Choice Pregnancy Midstream Test	Comments
Performance - Bench Top Used Sample Inv. ID# 242895	Positive Results Used HCG Hormone Positive Control 20miu/mL Lot 060710P Exp 6/20/10 (supplied by SCI)	Using a 5 second dip "urine" entered result window at 56 seconds; first reaction at 75 seconds; complete results at 80 seconds; constant for 3 minutes	Using a 10 second dip "urine" entered round window at 23 seconds; first reaction at 28 seconds; "urine" entered square window; first reaction at 27 seconds complete positive result at 48 seconds; constant for 10 minutes	Used SCI's instructions
	Negative Results Used Quantimetrix The Dipper Lot: 44501 Exp 2011-01	Using a 5 second dip "urine" entered result window at 51 seconds; first reaction at 75seconds; complete results at 85 seconds; constant for 3 minutes	Using a 10 second dip "urine" entered round window at 21 seconds; first reaction at 26 seconds; "urine" entered square window; first reaction at 31 seconds; complete negative results at 40 seconds; constant for 10 minutes	Used SCI's instructions
Physical	DIMENSIONS, inches Overall Length Width at Cap Width at Handle	5-5/8 13/16 3/4	6 11/16 7/16 narrowest	
	ABSORBENT TIP Length Width Thickness	1-1/16 7/16 1/8	1-1/2 (at tapered end of stick) 1/2 1/8	
	HANDLE Length Width Thickness	4-3/8 13/16 5/16	4-1/2 (tapered) 11/16 - 7/16(center) - 9/16 5/16	
	CAP Length Width Thickness	1-5/8 1-3/16 5/16	2-1/4 - 1-7/8 (tapered) 5/8 - 3/4 slightly tapered -5/8 5/16	

		First Response Pregnancy Test Early Result (6 days before a missed period)	SCI Brand New Choice Pregnancy Midstream Test	Comments
Physical	Construction	Opaque/white plastic contoured handle with clear/pink plastic tapered cap; white rectangular-shaped absorbent tip at one end. At mid section, one oval result window covered with plastic with printed results and symbols to one side of it; thumb grips at end of handle; First Response and symbol embossed on lower handle.	Opaque/white plastic tapered stick with pink tapered plastic cover; white rectangular-shaped absorbent tip at one end. At mid section of handle, one result window with 2 sections: one round, one irregular shape with printed results and symbols to one side of it; grip thumb line near tip end and several lines at end (samples submitted)	
	Product Workmanship	Acceptable	Acceptable	
	Packaging	Individually wrapped in printed heat-sealed foil pouch; boxed with clear plastic overwrap. Instruction pamphlet with illustrations inside.	Individually wrapped sticks in heat-sealed foil packets (use instructions submitted)	