



Specialized Technology Resources
 85 John Road
 Canton, MA 02021
 Phone: 781-821-2200
 Fax: 781-821-9266
 www.STRQuality.com

Technically Advanced Quality Assurance TAQA® Evaluation Report

Mr. Mahmoud M. Sales
SCI International, Inc.
 5902 Enterprise Court
 Fredrick, MD 21703

Folder #: 1013376-2 & 1015170
 Date: 6/25/2010 & 7/23/2010
 Received Date: 6/8/10 & 7/1/10
 Page: 1 of 3

cc: Mr. Stuart Neuman

| | | |
|-----------|---|---------------------------------------|
| Product: | SCI Brand New Choice Pregnancy Midstream Test* | UPC Code: |
| Size: | Sample | Reference #: Prepaid American Express |
| Supplier: | SCI International, Inc. | Inventory ID #: 239460 & 242895* |
| Code: | Lot 01209C03 Exp 04/2011 and Lot 0620C03 Exp 07/2012* | Lab Sample #: 1013376001 & 1015170001 |
| Purpose: | TAQA EVALUATION vs. E.P.T. EARLY RESPONSE | |

National Reference Brand: E.P.T. Early Pregnancy Test Result, (±) Results 5 days Before Missed Period;
 Code: Exp 2012 04 Lot No. 27092 - provided by SCI International, Inc.

Overall Rating: 10.0 Sample submitted meets TAQA quality criteria for parameters examined.

Note: * SCI Brands includes: New Choice, Perfect Choice, Baby Test, Baby Check, UCheck

Overall Quality Rating *

- 9.0 - 10.0** Sample tested is Comparable to the National Brand; meets TAQA quality criteria
- 7.5 - < 9.0** Sample tested is Marginally Comparable to the National Brand; does not meet all TAQA quality criteria
- < 7.5** Sample tested is Not Comparable to the National Brand; does not meet TAQA quality criteria

 Digitally signed by Maryann Sullivan
 Date: 2010.07.23 13:43:27 -04'00'

Maryann W. Sullivan, Client Services Manager

| Determination | | E.P.T. Early Pregnancy Test (±) Results 5 days Before Missed Period | SCI Brand New Choice Pregnancy Midstream Test | Comments |
|---|--|---|--|-----------------------------|
| Performance - Bench Top Used Sample Iv. # 242895 | Positive Results Used HCG Hormone Positive Control 20miu/mL Lot 060710P Exp 6/20/10 (supplied by SCI) | Using a 20 sec. dip "urine" entered the round window at 13 sec. with a first reaction at 28 sec.; "urine" entered square window at 29 sec. with results at 41 seconds; complete positive at 53 sec.; constant for 10 minutes | Using a 10 second dip "urine" entered round window at 23 seconds, first reaction at 28 seconds; "urine" entered square window, first reaction at 27 seconds; complete positive result at 48 seconds; constant for 10 minutes | Used SCI's use instructions |
| | Negative Results Used Quantimetrix The Dipper Lot: 44501 Exp 2011-01 | Using a 20 second dip "urine" entered round window at 20sec. with first reaction at 29sec.; "urine" entered square window at 29 sec. with a results at 44 sec.; complete negative results at 55ec.; constant for 10 minutes | Using a 10 second dip "urine" entered round window at 21 seconds, first reaction at 26 seconds, "urine" entered square window, first reaction at 31 seconds; complete negative results at 40 seconds; constant for 10 minutes | Used SCI's use instructions |
| Physical | DIMENSIONS, inches Overall Length Width at Cap Width at Handle | 5-1/2 1 7/8 | 6 11/16 7/16 narrowest | |
| | ABSORBENT TIP Length Width Thickness | 1 (rounded end) 1/2 1/8 | 1-1/2 (at tapered end of stick) 1/2 1/8 | |
| | HANDLE Length Width Thickness | 4-3/8 1 (widest) - 5/16 (tapered) 3/8 (widest) | 4-1/2 (tapered) 11/16 - 7/16 (tapered) 5/16 | |
| | CAP Length Width Thickness | 2-1/4 - 1-3/4 (tapered) 5/8 - 7/8 slightly tapered 5/16 | 2-1/4 - 1-7/8 (tapered) 5/8 - 3/4 slightly tapered -5/8 5/16 | |

| | | E.P.T. Early Pregnancy Test (±) Results 5 days Before Missed Period | SCI Brand New Choice Pregnancy Midstream Test | Comments |
|-----------------|---------------------|--|--|----------|
| Physical | Construction | Opaque/white plastic contoured handle with opaque/purple plastic tapered cap; white rectangular-shaped absorbent tip with rounded end. At mid section, one oval result window covered with plastic with 2 sections: one round, one rectangular; thumb grips line near tip end, several raised lines at end of handle; "e.p.t." printed in purple on handle | Opaque/white plastic tapered stick with pink tapered plastic cover; white rectangular-shaped absorbent tip at one end. At mid section of handle, one result window with 2 sections: one round, one irregular shape with printed results and symbols to one side of it; grip thumb line near tip end and several lines at end (samples submitted) | |
| | Product Workmanship | Acceptable | Acceptable | |
| | Packaging | Individually wrapped in printed heat-sealed foil pouch; boxed with clear plastic overwrap. Instruction pamphlet with illustrations inside. | Individually wrapped sticks in foil heat-sealed packets (use instructions submitted) | |